



GENERAL LINE **TRACTION CATALOG**



CONTENTS

| | |
|---------------------------------------|---|
| Introduction | 1 |
| Traction treatment: a background..... | 2 |
| Galaxy* TTET400 Traction Table | 3 |
| Galaxy* TTET300 Traction Table..... | 4 |
| Traction Kits | 5 |
| TX™ Traction Unit | 6 |
| Traction Accessories..... | 7 |

TRACTION TREATMENT: A BACKGROUND

In modern medicine, spinal distraction has been used as a treatment option to relieve back pain resulting from herniated, bulging, or protruding discs associated with various spine injuries and pathologies.

Stretching the spine helps reduce pressure on the intervertebral disc, relieves compression and irritation of the nerve roots and improves sagittal spine alignment.

CLINICAL INDICATIONS FOR TRACTION THERAPY

The traction device helps provide traction and mobilization of skeletal structures and skeletal muscles. The traction device may be used to help relieve peripheral radiation/sciatica and pain associated with:

- Protruding discs
- Bulging discs
- Herniated discs
- Degenerative disc disease
- Posterior facet syndrome
- Acute facet problems
- Radicular pain
- Prolapsed discs
- Spinal root impingement
- Hypomobility
- Degenerative joint disease
- Facet syndrome
- Compressions fracture
- Joint pain
- Discogenic pain



GALAXY*

TTET400 TRACTION TABLE



*table shown with optional axilla post package

GALAXY TTET400 4 SECTION VALUE, FORM AND FUNCTION

The Galaxy TTET400 is a multifunctional 4-section traction table which, in conjunction with a traction unit, aids in efficient spinal traction treatment. The TTET400's roller-mounted surface helps prevent friction below the hips and legs, aiding in full traction force during vertebral treatment. In addition, because the rolling sections can be locked in place, this variable-height traction table can also serve as an additional treatment table.

FEATURES

- A Scissor frame with four feet and four multidirectional, retractable wheels for easy transport
 - Provides maximum patient comfort in prone and supine positions to aid both manual therapy and traction therapy
 - Roller-mounted surface helps prevent friction below the hips and legs
 - Repositional head /back section elevates 82° and leg section to 73°
 - Gas spring assists on friction-free moveable sections for easy positioning
 - HalloTronic™+ System enables electric high/low height adjustment from 19 to 39 cm with capacity up to 440 LBs
- *Optional axilla post package for added stability for your patient. Must be added at the time of order.

WARRANTY

DJO Global® provides a 24-month warranty on the delivered product and a 12-month warranty on its upholstered parts, starting from the date of original consumer purchase, as indicated in a proof of purchase. See IFU for more detail.

* GALAXY is a table name specific to DJO® tables manufactured by MEDEN-INMED exclusively for DJO
+Hallotronic™ is a trademark by MEDEN-INMED Spółka z o.o.

ORDERING INFORMATION

| PART NO. | DESCRIPTION |
|----------|-------------|
|----------|-------------|

GALAXY TTET400

32671XUS Scissor frame with 4 feet and 4 multidirectional, retractable wheels

OPTIONAL AXILLA POST PACKAGE

37185 Axilla Posts Option is a special order and can be added at a later date. This option includes the following:
2 Axilla Posts, Yoke installed on upper frame, Docking elements installed on traction arm

(X) UPHOLSTERY COLOR OPTIONS

| | | |
|---------------------------|-------------------|---------------------------|
| Beige (x = 01) | Blue (x = 02) | Gray (x = 06) |
| Graphite Gray (x = 07) | Black (x = 08) | Imperial Blue (x = 19) |

When ordering, replace 'x' in Part Number with desired 2-digit color code above. For example, 3267101 = Galaxy TTET400, beige

NOTICE: Color charts are produced within the limitations of lithographic color printing and are intended as a general guide only.

TECHNICAL SPECIFICATIONS

| | | |
|---------------------------------|-------------------------------------|--------------------------------------|
| Mains Power: | 120 V~ 50/60 HZ | |
| Current: | Max 2A with 120V | |
| Dimensions: | Head section | 23 inches |
| | Middle section | 8 inches |
| | Pelvic section | 9.8 inches |
| | Leg section | 47 inches |
| | Width | 27.5 inches |
| | Total table length | 79.5 inches |
| Adjustment: | Total length with traction platform | 103 inches |
| | Angle Inclination | Head: -10° to +82° Leg: 0° to 73° |
| | Displacement | Leg: 6 inches |
| | Height | min: 19 inches max: 39 inches |
| Dynamic lifting capacity: | 440 LBs | |
| Max. workload: | 605 LBs | |
| Electrical Safety Class: | Class I, Type B | |
| Safety Tests: | CE, Class I Rule 12 | |
| Shipping Weight | 330 LBs | |
| Shipping Dimns. (LxWxH) inches: | 86.6x33.5x31.5 inches | |
| Table Weight: | 274 LBs | |

* Dimensions are listed as L x W x H

ORDERING INFORMATION

| PART NO. | DESCRIPTION |
|----------|-------------|
|----------|-------------|

GALAXY TTET300

| | |
|----------|--|
| 32641xUS | Scissor frame with 4 feet and 4 multidirectional, retractable wheels |
|----------|--|

OPTIONAL AXILLA POST PACKAGE

| | |
|-------|--|
| 37185 | Axilla Posts Option is a special order and can be added at a later date. This option includes the following: 2 Axilla Posts, Yoke installed on upper frame, Docking elements installed on traction arm |
|-------|--|

(X) UPHOLSTERY COLOR OPTIONS

| | | |
|---------------------------|-------------------|---------------------------|
| Beige (x = 01) | Blue (x = 02) | Gray (x = 06) |
| Graphite Gray (x = 07) | Black (x = 08) | Imperial Blue (x = 19) |

When ordering, replace 'x' in Part Number with desired 2-digit color code above. For example, 3264101 = Galaxy TTET300, beige

NOTICE: Color charts are produced within the limitations of lithographic color printing and are intended as a general guide only.

TECHNICAL SPECIFICATIONS

| | | |
|---------------------------------|-------------------------------------|----------------------------------|
| Mains Power: | 120 V~ 50/60 HZ | |
| Current: | Max 2A with 120V | |
| Dimensions: | Head section | 23 inches |
| | Middle section | 8 inches |
| | Leg section | 44.7 inches |
| | Width | 27.5 inches |
| | Total table length | 79.5 inches |
| | Total length with traction platform | 103 inches |
| Adjustment: | Angle Inclination | Head: -10° to +82° |
| | Displacement | Leg: 6 inches |
| | Height | min: 19 inches max: 39 inches |
| Dynamic lifting capacity: | 440 LBs | |
| Max. workload: | 605 LBs | |
| Electrical Safety Class: | Class I, Type B | |
| Safety Tests: | CE, Class I Rule 12 | |
| Shipping Weight | 330 LBs | |
| Shipping Dimns. (LxWxH) inches: | 86.6x33.5x31.5 inches | |
| Table Weight: | 274 LBs | |

* Dimensions are listed as L x W x H

GALAXY*

TTET300 TRACTION TABLE



*table shown with optional axilla post package

GALAXY TTET300 3 SECTION VALUE, FORM AND FUNCTION

The Galaxy TTET300 is a multifunctional 3-section traction table which, in conjunction with a traction unit, aids in efficient spinal traction treatment. The TTET300's roller-mounted surface helps prevent friction below the hips and legs, aiding in the full traction force during vertebral treatment. In addition, because the rolling sections can be locked in place, this variable-height traction table can also serve as an additional treatment table.

FEATURES

- A scissor frame with four feet and four multidirectional, retractable wheels for easy transport
- Provides maximum patient comfort in prone and supine positions to aid both manual therapy and traction therapy
- Roller-mounted surface helps prevent friction below the hips and legs
- Repositional head /back section elevates to 82°
- Gas spring assists on friction-free moveable sections for easy positioning
- HalloTronic™ + System enables electric high/low height adjustment from 19 to 39 inches with capacity up to 440 LBs

*Optional axilla post package for added stability for your patient. Must be added at the time of order.

WARRANTY

DJO Global® provides a 24-month warranty on the delivered product and a 12-month warranty on its upholstered parts, starting from the date of original consumer purchase, as indicated in a proof of purchase. See IFU for more detail.

* GALAXY is a table name specific to DJO® tables manufactured by MEDEN-INMED exclusively for DJO

+Hallotronic™ is a trademark by MEDEN-INMED Spółka z o.o

TX™ TRACTION UNIT

FEATURES

TX™ Traction Unit

- Static, intermittent, and cycling traction
- User-defined hold, rest and treatment times
- 10 User-defined protocol slots to store customized settings
- Digital monochromatic touch screen interface that pivots 270°
- Easily store treatment parameters, pain scales, and pain maps on Patient Data Cards



TRACTION KITS



ORDERING INFORMATION

| PART NO. | DESCRIPTION |
|--------------------------------------|---|
| TRU-TRAC™ TRACTION DEVICE | |
| 4759 | TX™ Traction Unit (including patient switch and data card) |
| TECHNICAL SPECIFICATIONS | |
| Mains Power: | 100 - 240 V, 50/60 Hz |
| Power Consumption: | 75 W |
| Current Consumption: | 3.2 Amps Max |
| Weight: | 30 lbs |
| Dimensions*: | 18 x 9 x 18 inches |
| Electrical Safety Class: | Class 1, Type B |
| Safety Tests: | IEC/UL/EN 60601-1, IEC/EN 60601-1-2, CAN C22.2 No. 601.1-M90 w/A2 - Meets MDD 93/42/EEC, CE 0413 |
| * Dimensions are listed as L x W x H | |

ORDERING INFORMATION

| PART NO. | DESCRIPTION |
|---|---|
| TRACTION KITS WITH STOOLS | |
| 8200 | Traction Kit - Beige Stool |
| 8201 | Traction Kit - Blue Stool |
| 8202 | Traction Kit - Grey Stool |
| 8203 | Traction Kit - Graphite Grey Stool |
| 8204 | Traction Kit - Black Stool |
| 8205 | Traction Kit - Imperial Blue Stool |
| TRACTION KITS LISTED ABOVE INCLUDE THE FOLLOWING STANDARD ACCESSORIES: | |
| 4759 | TX™ Traction Unit |
| A. 2825 | QuikWrap Belting System |
| C. 3265* | Traction Stool (6 color choices) |
| D. 7040 | Saunders Cervical Traction System - includes Clevis |
| 8206 - TX™ TRACTION KIT WITH BOLSTERS | |
| 4759 | TX™ Traction Unit |
| A. 2825 | QuikWrap Belting System |
| B. 42901 | Knee Bolster |
| D. 7040 | Saunders Cervical Traction System - includes Clevis |

ORDERING INFORMATION

| PART NO. | DESCRIPTION |
|---|---------------------------------|
| A QUIKWRAP BELTING SYSTEM | |
| 2825 | QuikWrap Belt Set |
| B KNEE BOLSTER | |
| B1 42901 | Knee Bolster |
| B2 42904 | Knee Bolster |
| C KALTENBORN | |
| 371308 | Kaltenborn 3.7x6x2 inches |
| D AXIAL POSTS | |
| D 83215ME | Traction Axial Post, Black (x2) |
| <p>The following items are available in 6 colors: Beige (01) Blue (02) Grey (06) Graphite Grey (07) Black (08) Imperial Blue (19). Reference * for the 6 color choices below.</p> | |
| E WEDGE | |
| E1 3714* | 11.8x15.7x4.7 inches |
| E2 3715* | 19.7x13.7x7.9 inches |
| F ROUND BOLSTER | |
| 3718* | Round Bolster 23.6x6 inches |


(X) WEDGE AND BOLSTER COLOR OPTIONS

| | | |
|---------------------------|-------------------|---------------------------|
| Beige (x = 01) | Blue (x = 02) | Gray (x = 06) |
| Graphite Gray (x = 07) | Black (x = 08) | Imperial Blue (x = 19) |

When ordering, replace 'x' in Part Number with desired 2-digit colour code above. For example, 371808 = Black Round Bolster

NOTICE: Color charts are produced within the limitations of lithographic color printing and are intended as a general guide only.

TRACTION ACCESSORIES



FEATURES

A QuikWrap Belt Set

- Universal fit, for use with pelvic or lumbar traction applications up to 300 LBs

B Knee Bolster

- Provides optimum body position to aid in multiple treatment applications.

C Kaltenborn Wedge

- For use by professionals trained in manual mobilization techniques.

D Axial Posts

- Removable axilla posts help provide additional fixation to the table, preventing slippage during supine traction

NOTE: These posts will only work with the Galaxy TTET traction tables if the Axilla Post Package has been added at the time of purchase with the table.

E Wedge (1 per pack)

- Flat surfaced positioning shapes in two sizes 11.8x15.7x4.7 or 19.7x13.7x7.9 inches

F Round Bolster

- Rounded positioning shape in one size: 23.6x6 inches



G Facial Pillow

- Supports the head and neck of the patient for more comfortable traction sessions

H Pneumatic Stools

- High-quality, therapy stools with 60 Kg/m³ foam density swivel seat, supports 330 LBs

I Traction Flexion Stool

- Adjustable height enabling placement of the patient's lumbar spine into a neutral or flexed position for increased comfort. Upholstery can be color coordinated to match traction table

J Pron Pillo®

- Elevates the head and shoulders for added comfort. The raised base helps allow patients in prone position to breathe more easily

K Stabilizer Pressure Biofeedback

- Consists of combined gauge and inflation bulb connected to pressure cell

L Clinical Cervical Traction

- This clinical device is designed specifically to aid in cervical traction. One size fits most

M Safety Key

- Removable personal safety key for table security and positioning.

ORDERING INFORMATION

| PART NO. | DESCRIPTION |
|---|---|
| G FACIAL PILLOW | |
| 4291 | Facial Pillow |
| H PNEUMATIC STOOLS (See color chart on page 7) | |
| 371901 | Standard Stool Beige |
| 371902 | Standard Stool Blue |
| 371906 | Standard Stool Gray |
| 371907 | Standard Stool Graphite Gray |
| 371908 | Standard Stool Black |
| 371919 | Standard Stool Imperial Blue |
| 372001 | Flexi Stool Beige |
| 372002 | Flexi Stool Blue |
| 372006 | Flexi Stool Gray |
| 372007 | Flexi Stool Graphite Gray |
| 372008 | Flexi Stool Black |
| 372019 | Flexi Stool Imperial Blue |
| I TRACTION FLEXION STOOLS (See color chart on page 7) | |
| 3265101 | Beige |
| 3265102 | Blue |
| 3265106 | Gray |
| 3265107 | Graphite Gray |
| 3265108 | Black |
| 3265119 | Imperial Blue |
| J PRON PILLO | |
| 4206 | Pron Pillo |
| K STABILIZER PBF | |
| 9296 | Stabilizer Pressure Biofeedback (PBF) Unit |
| STANDARD ACCESSORIES | |
| 92965 | Instruction Manual |
| 92961 | PBF Stabilizer Pressure Cuff |
| L SAUNDERS CERVICAL TRACTION | |
| 7040 | Saunders Cervical Traction System - Includes Clevis |
| M SAFETY KEY | |
| 15-0127 | Safety Key |



www.chattanoogaarehab.com

Triton® and Galaxy Tables are manufactured by Meden-Inmed Sp. z o.o. exclusively for Chattanooga. A DJO Global Company